# TD.01533

#### Rodent Diet (2018, 6000 Doxycycline)

Formula	g/Kg
2018, Teklad Global 18% Protein Rodent Diet	994.0
Doxycycline Hyclate	6.0
	2018, Teklad Global 18% Protein Rodent Diet



# **Key Features**

- Standard Diet Base
- 2018
- Doxycycline
- Gene Induction

# **Key Planning Information**

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week More frequent replacement may be advised
- Lead time:
  - · 2 weeks non-irradiated
  - · 4 weeks irradiated



# **Footnote**

Calculated for 6000 mg/kg doxycycline hyclate (approx. 5200 ppm doxycycline, as the hyclate form contains ~87% doxycycline).

# **Product Specific Information**

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

# Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## Place Your Order (USA & Canada)

Please Choose One

% kcal from

www.inotivco.com/teklad-orders

**Price Quotes | Order Status** 

- tekladorders@inotivco.com
- 800.483.5523

**Selected Nutrient Information** 

70 Dy Weigiit	/o KCai II OIII
18.5	24.4
43.9	57.9
6.0	17.7
	18.5 43.9

3.0 Values are calculated from ingredient analysis or manufacturer data

<sup>2</sup>Estimated digestible carbohydrate

Kcal/g

#### Collaborate with a Nutritionist

askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

teklad@inotivco.com