

# TD.03142

## 8% NaCl Diet (Red)

Formula	g/Kg
Wheat	350.0
Corn	214.49
Soybean Meal	201.0
Corn Gluten Meal (60%)	56.0
Alfalfa Meal (17%), dehydrated	30.0
Corn Oil	35.8
Sodium Chloride	80.0
Dicalcium Phosphate, FG (18.5% P, 21% Ca)	14.0
Calcium Carbonate, FG (38%)	12.0
Mineral Mix, TSD (80318)	1.5
Vitamin Mix, TSD (81125)	3.0
DL-Methionine, FG (99%)	1.0
L-Lysine HCl, FG (78%)	1.0
Ethoxyquin, antioxidant	0.01
Red Food Color	0.2

### Footnote

This formula is a modification of TD.92012 to add red food coloring. This formula is part of a series of adjusted NaCl formulas of similar composition. Approximate nutrient content: 19% protein, 5% fat, 0.9% Ca, 0.65% P, 0.15% Mg, 0.8% K, 3.1% Na, and 5% Cl.

### Selected Nutrient Information<sup>1</sup>

	% by weight	% kcal from
Protein	18.5	25.3
Carbohydrate <sup>2</sup>	42.8	58.4
Fat	5.3	16.4
Kcal/g	2.9	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

<sup>2</sup> Estimated digestible carbohydrate

### Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)



### Key Features

- Natural Ingredient Diet
- Adjusted NaCl Series
- Color Coded
- Rodent

### Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



### Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

### Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

### Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)