

# TD.03314

## 62% Sucrose Diet (No Fat)

Formula	g/Kg
Casein, "Vitamin-Free" Test	219.0
DL-Methionine	2.5
Sucrose	620.0
Cellulose	109.5
Vitamin Mix, Teklad (40060)	10.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0

### Footnote

This is a modification of TD.93165 to remove the corn oil (5%) and increase sucrose (2%) and cellulose (3%). Pellets are hard. Diet is free of total fat and essential fatty acids so should be used for specific research purposes only.

### Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	20.0	24.1
Carbohydrate	62.8	75.9
Fat	0.0	0.0
Kcal/g	3.3	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

### Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only  
Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



### Key Features

- Purified Diet
- Sucrose
- Fat Free
- Rodent

### Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



### Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

### Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

### Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)