<b>TD.05298</b> Dox Diet (2018, 1000)			inotiv
			analyze. answer. advance
Formula		g/Kg	Key Features
2018, Teklad Global 18% Protein Doxycycline Hyclate	Rodent Diet	999.0 1.0	<ul> <li>Standard Diet Base</li> <li>Doxycycline</li> <li>Gene Induction</li> <li>Teklad 2018</li> </ul>
			<ul> <li>Key Planning Information</li> <li>Products are made fresh to order</li> <li>Store product at 4°C or lower</li> <li>Use within 6 months (applicable to most dialed)</li> <li>Box labeled with product name, manufacturing date, and lot number</li> <li>Replace diet at minimum once per week More frequent replacement may be advise</li> <li>Lead time: <ul> <li>2 weeks non-irradiated</li> <li>4 weeks irradiated</li> </ul> </li> </ul>
Footnote			Product Specific Information
Calculated for 1000 mg/kg doxycycline hyclate (approx. 860 ppm doxycycline, as the hyclate form contains ~86% doxycycline).			<ul> <li>1/2" Pellet or Powder (free flowing)</li> <li>Minimum order 3 Kg</li> <li>Irradiation available upon request</li> </ul>
			Options (fees will apply) <ul> <li>Rush order (pending availability)</li> <li>Irradiation (see Product Specific Information</li> <li>Vacuum packaging (1 and 2 Kg)</li> </ul>
Selected Nutrient Infor	mation <sup>1</sup>		Place Your Order (USA & Canad
<b>.</b>	% by weight	% kcal from	Please Choose One
Protein Carbohydrate <sup>2</sup>	18.6 44.2	24.4 57.9	<ul> <li>www.inotivco.com/teklad-orders</li> <li>tekladorders@inotivco.com</li> </ul>
Fat	6.0	17.7	• 800.483.5523
Kcal/g 3.0			
<sup>1</sup> Values are calculated from ingredient	analysis or manufacturer data		
<sup>2</sup> Estimated digestible carbohydrate			
Collaborate with a Nutritionist			Price Quotes   Order Status
<ul> <li>askanutritionist@inotivco.c</li> </ul>	om		teklad@inotivco.com

Teklad<sup>™</sup> Custom Diets, Madison, WI | inotivco.com