

TD.07011

29% Fat Diet

Formula

Formula	g/Kg
Casein	289.0
DL-Methionine	3.33
Maltodextrin	207.3
Sucrose	90.5
Hydrogenated Vegetable Oil	274.1
Corn Oil	16.0
Cellulose	53.12
Vitamin Mix, Teklad (40060)	13.33
Mineral Mix, AIN-76 (170915)	46.66
Calcium Carbonate	6.66

Footnote

Modification of TD.93075 to replace corn starch with maltodextrin for pelleting.
Hydrogenated vegetable oil is added as a source of dietary trans fats.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	25.5	21.0
Carbohydrate	29.9	24.6
Fat	29.3	54.4
Kcal/g	4.9	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only
Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Diet Induced Obesity
- Trans Fatty Acids

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com