

# TD.07573

Dox Diet (2014, 625, Y)

## Formula

	g/Kg
2014, Teklad Global 14% Protein Rodent Maintenance Diet	998.975
Doxycycline Hyclate	0.625
Yellow Food Color	0.4

## Footnote

The 2014 diet with 625 mg doxycycline hyclate per kg diet, and color coded yellow. Doxycycline hyclate contains approximately 87% doxycycline. A mouse consuming 4 g/day would receive approx. 2.2 mg of doxycycline.

## Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	14.2	19.3
Carbohydrate <sup>2</sup>	51.4	70.2
Fat	3.4	10.4

**Kcal/g 2.9**

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

<sup>2</sup> Estimated digestible carbohydrate

## Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



## Key Features

- Standard Diet Base
- Gene Induction
- Doxycycline
- Tekald 2014

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



## Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

## Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)