

TD.09651

Doxycycline Diet (625 mg/kg, 2019)

Formula

g/Kg

2019, Teklad Global 19% Protein Rodent Diet, pre-extruded
Doxycycline Hyclate

999.375

0.625

Footnote

This diet was designed to deliver a daily dose of 2-3 mg of doxycycline based on consumption of 4-5 g/d by a mouse. Doxycycline hyclate contains approximately 87% doxycycline. The diet is not extruded.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	18.8	21.4
Carbohydrate ²	51.5	58.7
Fat	7.8	19.9
Kcal/g	3.5	

¹ Values are calculated from ingredient analysis or manufacturer data

² Estimated digestible carbohydrate

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com

inotiv

analyze. answer. advance.

Key Features

- Standard Diet Base
- Doxycycline
- Gene Induction
- 2019

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com