

# TD.120658

Rodent Diet (2018, 1000 Dox, B)

## Formula

|   | g/Kg  |
|---|-------|
| 2018, Teklad Global 18% Protein Rodent Diet | 998.6 |
| Doxycycline Hyclate                         | 1.0   |
| Blue Food Color                             | 0.4   |

## Footnote

Modification of TD.05298 to add blue coloring. This diet was designed to deliver a daily dose of 3.5-4.4 mg of doxycycline based on consumption of 4-5 g/d by a mouse. Doxycycline hyclate contains approximately 87% doxycycline.

## Selected Nutrient Information <sup>1</sup>

|                           | % by weight | % kcal from |
|---------------------------|-------------|-------------|
| Protein                   | 18.6        | 24.4        |
| Carbohydrate <sup>2</sup> | 44.1        | 57.9        |
| Fat                       | 6.0         | 17.7        |
| <b>Kcal/g</b>             | <b>3.0</b>  |             |

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

<sup>2</sup> Estimated digestible carbohydrate

## Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



## Key Features

- Standard Diet Base
- 2018
- Doxycycline
- Color Coded Blue

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



## Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

## Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)