

TD.120769

Rodent Diet (2018, 625 Dox, B)

Formula

	g/Kg
2018, Teklad Global 18% Protein Rodent Diet	998.975
Doxycycline Hyclate	0.625
Blue Food Color	0.4

Footnote

Modification of TD.01306 to add blue coloring. This diet was designed to deliver a daily dose of 2-3 mg of doxycycline based on consumption of 4-5 g/d by a mouse. Doxycycline hyclate contains approximately 87% doxycycline.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	18.6	24.4
Carbohydrate ²	44.2	57.9
Fat	6.0	17.7
Kcal/g	3.0	

¹ Values are calculated from ingredient analysis or manufacturer data

² Estimated digestible carbohydrate

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Standard Diet Base
- 2018
- Doxycycline
- Color Coded Blue

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com