

TD.130141

Doxycycline Diet (2020, 625 Dox, R)

Formula	g/Kg
2020, Teklad Global Soy Protein-Free Rodent Diet, pre-extruded	999.175
Doxycycline Hyclate	0.625
Red Food Color	0.2

Footnote

Modification of TD.110720 to add red coloring. This diet was designed to deliver a daily dose of 2-3 mg of doxycycline based on consumption of 4-5 g/d by a mouse. Doxycycline hyclate contains approximately 87% doxycycline. The diet is not extruded.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	19.0	23.7
Carbohydrate ²	48.9	61.2
Fat	5.4	15.1

Kcal/g 3.2

¹ Values are calculated from ingredient analysis or manufacturer data

² Estimated digestible carbohydrate

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Standard Diet Base
- 2020
- Doxycycline
- Color Coded Red

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com