TD.180327

40 ppm Doxy Diet (2018, Red)

Formula	g/Kg
2018, Teklad Global 18% Protein Rodent Diet	999.704
Doxycycline Hyclate	0.046
Red Food Color	0.25



Key Features

- Standard Diet Base
- 2018
- Doxycycline
- Red Food Color

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week More frequent replacement may be advised
- Lead time:
 - · 2 weeks non-irradiated
 - · 4 weeks irradiated



Footnote

Designed to add 40 mg doxycyline/kg diet using the Teklad 2018 Global Rodent Diet as the base. Doxycycline hyclate contains approximately 87% doxycycline. A red food dye is added.

Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Selected Nutrient Information ¹

		% by weight	% kcal from
Protein		18.6	24.4
Carbohydrate ²		44.2	57.9
Fat		6.0	17.7
Kcal/g	3.1		

3.1 Values are calculated from ingredient analysis or manufacturer data

²Estimated digestible carbohydrate

Collaborate with a Nutritionist

askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only Teklad™ Custom Diets, Madison, WI inotivco.com

Price Quotes | Order Status

teklad@inotivco.com