

TD.88232

Basal Mix (without 20% Oil)

Formula	g/Kg
Casein	300.0
DL-Methionine	4.5
Sucrose	373.5
Corn Starch	187.5
Cellulose	62.5
Mineral Mix, AIN-76 (170915)	52.5
Calcium Carbonate	4.5
Vitamin Mix, Teklad (40060)	15.0

Footnote

This basal mix is designed to be used at the rate of 800 g/Kg of diet, in conjunction with 200 g of the oil source(s).

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	26.6	32.1
Carbohydrate	55.5	67.1
Fat	0.3	0.8
Kcal/g	3.3	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only
Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Basal Mix
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com