TD.92163

Selenium Deficient Diet

| Formula | g/Kg |
|-----------------------------------|-------|
| Torula Yeast | 300.0 |
| DL-Methionine | 3.0 |
| Sucrose | 591.0 |
| Corn Oil | 50.0 |
| Mineral Mix, Se Deficient (80313) | 35.0 |
| Calcium Carbonate | 11.0 |
| Vitamin Mix, Teklad (40060) | 10.0 |
| 1 | |



Key Features

- Purified Diet
- Selenium
- Torula Yeast
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
 More frequent replacement may be advised
- Lead time:
 - · 2 weeks non-irradiated
 - · 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Selected Nutrient Information 1

| | % by weight | % kcal from |
|--------------|-------------|-------------|
| Protein | 15.3 | 16.5 |
| Carbohydrate | 64.1 | 69.2 |
| Fat | 5.9 | 14.3 |

Estimated background selenium for this kind of formula is less than 0.01 ppm Se. This is

slightly lower than for casein type diets. TD.96363 is a possible ingredient matched control

with selenium added at 0.1 mg/kg diet from sodium selenite. For questions about this or

Kcal/g 3.7

Footnote

Values are calculated from ingredient analysis or manufacturer data

other selenium adjusted diets askanutritionist@inotivco.com.

Collaborate with a Nutritionist

askanutritionist@inotivco.com

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

teklad@inotivco.com