## TD.93075

### Adjusted Fat Diet (55/Fat)

Formula	g/Kg
Casein	289.0
DL-Methionine	3.33
Corn Starch	207.3
Sucrose	90.5
Hydrogenated Vegetable Oil	274.1
Corn Oil	16.0
Cellulose	53.12
Vitamin Mix, Teklad (40060)	13.33
Mineral Mix, AIN-76 (170915)	46.66
Calcium Carbonate	6.66



## **Key Features**

- Purified Diet
- Diet Induced Obesity
- Trans Fatty Acids
- Dough Form

## **Key Planning Information**

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
  More frequent replacement may be advised
- Lead time:
  - · 2 weeks non-irradiated
  - · 4 weeks irradiated



# **Product Specific Information**

- Powder (dough)
- Minimum order 3 Kg
- Irradiation available upon request

#### Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)

# Selected Nutrient Information <sup>1</sup>

trans fats. Cited in Diabetes (2004) 53:3274-3285.

Science italicite illioilliation		
	% by weight	% kcal from
Protein	25.5	21.2
Carbohydrate	28.8	24.0
Fat	29.3	54.8

A diet used to induce obesity. Hydrogenated vegetable oil is added as a source of dietary

Kcal/g 4.8

**Footnote** 

<sup>1</sup>Values are calculated from ingredient analysis or manufacturer data

#### **Collaborate with a Nutritionist**

askanutritionist@inotivco.com

#### Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

#### Price Quotes | Order Status

• teklad@inotivco.com