

TD.94047

Vitamin Mix, AIN-93-VX

Formula	g/Kg
Niacin	3.0
Calcium Pantothenate	1.6
Pyridoxine HCl	0.7
Thiamin (81%)	0.6
Riboflavin	0.6
Folic Acid	0.2
Biotin	0.02
Vitamin B ₁₂ (0.1% in mannitol)	2.5
Vitamin E, DL-alpha tocopheryl acetate (500 IU/g)	15.0
Vitamin A Palmitate (500,000 IU/g)	0.8
Vitamin D ₃ , cholecalciferol (500,000 IU/g)	0.2
Vitamin K ₁ , phylloquinone	0.075
Sucrose, fine ground	974.705

Footnote

Reference: J. Nutr. 123:1939-1951, 1993. Supplies the recommended concentrations of vitamins for AIN-93G and AIN-93M diets, when used at 10 g/kg of diet. The choline source should be listed as a separate item in diets utilizing this vitamin mix.

Collaborate with a Nutritionist

- askanutritionist@inotivco.com



Key Features

- Vitamin Mix
- Commonly Used

Key Planning Information

- Products are made fresh to order
- Lead time 2 weeks
- Store product at 4°C or lower
- Use within 6 months (applicable to most mixes)
- Box labeled with product name, manufacturing date, and lot number



Product Specific Information

- Powder
- Minimum order 500 grams

Options (fees will apply)

- Rush order (pending availability)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com