

TD.94048

AIN-93M Purified Diet

Formula

g/Kg

Casein	140.0
L-Cystine	1.8
Corn Starch	465.692
Maltodextrin	155.0
Sucrose	100.0
Soybean Oil	40.0
Cellulose	50.0
Mineral Mix, AIN-93M-MX (94049)	35.0
Vitamin Mix, AIN-93-VX (94047)	10.0
Choline Bitartrate	2.5
TBHQ, antioxidant	0.008

Footnote

Reference: J. Nutr. 123:1939-1951, 1993. Formulated for the maintenance of adult, non-reproducing rodents.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	12.4	13.7
Carbohydrate	68.3	75.9
Fat	4.1	10.3
Kcal/g	3.6	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- American Institute of Nutrition
- Maintenance
- Commonly Used

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation not advised
Contact a nutritionist for recommendations

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com