

TD.95078

2% NaCl Diet



Formula	g/Kg
Wheat	350.0
Corn	289.49
Soybean Meal	194.1
Corn Gluten Meal (60%)	50.0
Alfalfa Meal (17%), dehydrated	30.0
Corn Oil	33.9
Sodium Chloride	20.0
Dicalcium Phosphate, FG (18.5% P, 21% Ca)	14.0
Calcium Carbonate, FG (38%)	12.0
Mineral Mix, TSD (80318)	1.5
Vitamin Mix, TSD (81125)	3.0
DL-Methionine, FG (99%)	1.0
L-Lysine HCl, FG (78%)	1.0
Ethoxyquin, antioxidant	0.01

Key Features

- Natural Ingredient Diet
- Adjusted NaCl Series
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com

Footnote

This diet is a modification of TD.90228. This formula is part of a series of adjusted NaCl formulas of similar composition. Approximate nutrient content: 19% protein, 5% fat, 0.9% Ca, 0.65% P, 0.15% Mg, 0.8% K, 0.8% Na, and 1.3% Cl. Contact Technical Services for details.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	18.4	23.6
Carbohydrate ²	47.5	60.9
Fat	5.4	15.5
Kcal/g	3.1	

¹ Values are calculated from ingredient analysis or manufacturer data

² Estimated digestible carbohydrate

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com