

TD.96121

21% MF, 1.25% Chol. Diet

Formula

g/Kg

Casein	195.0
DL-Methionine	3.0
Sucrose	341.46
Corn Starch	150.0
Anhydrous Milkfat	210.0
Cholesterol	12.5
Cellulose	39.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.04

Footnote

This is a modification of TD 88137 (21% anhydrous milkfat by weight) to increase the added cholesterol to 1.25%. Total cholesterol (added + from fat) is approximately 1.3%. Related modification that includes 0.5% cholic acid is TD.02028. For questions about this diet or modifications askanutritionist@inotivco.com

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	17.3	15.2
Carbohydrate	48.5	42.7
Fat	21.2	42.0
Kcal/g	4.5	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only
Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Atherogenic
- Cholesterol
- High Fat

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com